

LIGHT



Spirobank II™ Basic

Essential,
easy to use Spirometer

Ideal for **family doctors**,
occupational medicine, screening

Real-time tests via **USB**

Rechargeable, **long-life battery**

Traffic lights for immediate test interpretation

Backlit, **high-resolution display**

On screen results and curves preview



Always included: Winspiro PRO™ PC software with free update



MIR SpiroPRO™ New Spirometry software for **Windows and MacOS**



Available with both **DISPOSABLE**
or **REUSABLE** turbine flowmeter



Spirometry test: FVC, VC, IVC, IC ERV,
PRE/POST Bronchodilator comparison
with a wide range of selectable parameters

Minispir Light™

Simplified Computer-based
Spirometer

Real time COPD and Asthma screening
have never been so intuitive and inexpensive!

Spirometry test interpretation

Embedded Temperature sensor for BTPS
conversion



Always included: Winspiro Light™ simplified
PC software with free update



Available with **DISPOSABLE** turbine flowmeter



Spirometry test: FVC, VC, IVC, PRE/POST
Bronchodilator comparison