

MIR Spiro™ android iOS

Simple and intuitive App for Spirometry and Oximetry test directly on your Tablets



*A smarter way to perform spirometry!
Maximum ease of use in primary care*

An expert beside you

Highly innovative **Virtual Assistant** supporting the operator before, during and after the spirometry test with four main functions:


- ✓ didactic guide with recommendations for the correct execution of the test
- ✓ identification of the cause of errors during test execution
- ✓ suggestions both for the operator and for the patient to make an acceptable test
- ✓ automatic test interpretation

Sharing is caring

Changing our software experience, direct and simple as a tap, test can be printed, exported in pdf and sent by email.

The recorded data can be transferred in HL7 format to a web service for archiving tasks or for second opinions and are enabled for EMR access

MIR Spiro is an App for iPad, iPad mini, Android Tablets



Virtual Assistant

