

PERSONAL SMART

SMART ONE™



Smart One™ 



SmartOne

Peak Flow and FEV1  
now available on your Smart

Suitable from 5 to 93 years!

Ideal in the **self management** of Asthma, COPD, Lung Transplant care, Cystic Fibrosis and for use in Clinical Trials

Share results in pdf (WhatsApp, email etc.) and print directly.

**Software Development Kit (SDK) for App Developers**  
Implemented for iOS and Android.

Source code of sample app provided



Spirometry test: PEF, FEV1




Available with **REUSABLE**  
single patient turbine flowmeter

Unparalleled Accuracy Unmatched Quality!


ISO 26782 Std (for spirometry parameters)

ISO 23747 Std (only for the Peak Flow parameter)

ATS/ERS 2019 standards



Traffic light health indicator for PEF, symptoms scoring and notes



Graphic trend of Peak Flow or FEV1