

ADVANCED



Spirobank II™ Advanced

Spirometer with Oximetry Option

Portable minilab for **pulmonologists and respiratory therapists**

Real-time tests with **wireless Bluetooth™** and **USB**

Rechargeable, **long-life battery**

Traffic lights for immediate test interpretation

SpO2 and pulse rate (optional) directly on the display (including **plethysmographic curve**)



Always included: Winspiro PRO™ PC software with free update



MIR SpiroPRO™ New Spirometry software for **Windows and MacOs**



Available with both **DISPOSABLE** or **REUSABLE** turbine flowmeter



Spirometry test: FVC, VC, IVC, MVV, PRE/POST Bronchodilator comparison with a wide range of selectable parameters



Optional function available:
Oximeter with adult or paediatric finger probe

Minispir™

Computer-based Spirometer

Real time Flow/Volume and Volume/Time curves on your PC for a comprehensive Spirometry

Includes a wide range of selectable parameters

Spirometry test interpretation

Embedded Temperature sensor for BT/PS conversion



Always included: Winspiro PRO™ PC software with free update



MIR SpiroPRO™ New Spirometry software for **Windows and MacOs**



Available with both **DISPOSABLE** or **REUSABLE** turbine flowmeter



Spirometry test: FVC, VC, IVC, MVV, PRE/POST Bronchodilator comparison with a wide range of selectable parameters